

iTero Restorative Material Selection Cheat Sheet

Ragle Product	iTero Universal Generic Restorative Material
Cera-Leux™ HSZ	Ceramic: Full Contour Zirconia **
Cera-Leux™ HSZ PLUS	Ceramic: Full Contour Zirconia **
Cera-Leux™ HTZ	Ceramic: Full Contour Zirconia **
Cera-Leux™ HTZ PLUS	Ceramic: Full Contour Zirconia **
Cera-Leux™ e.max®	Ceramic: Lithium Disilicate Milled**
Cera-Leux™ e.max® PLUS	Ceramic: Lithium Disilicate Milled**
Milled Full Gold High Noble (yellow)	Metal/PFM: Full Contour Gold **
Milled Full Gold Noble (yellow)	Metal/PFM: Full Contour Gold **

Step 1: Select your iTero Universal Definition Material based on the above chart.

** **Step 2:** Verify the specific Ragle Product Brand and/or alloy you are requesting by typing in the Notes section your product/alloy (i.e. “HSZ PLUS”, “FGC Noble”, etc.).

Step 3: If you would like the product’s “PLUS” option, please write this in the Notes. The PLUS options are esthetically enhanced with micro hand-layering, great for the anterior section of the mouth. *If not specified in the Notes, the PLUS option is the lab default for anterior teeth.*

